


## FT 09/1 Thermometer blue


 german  
 engineering

**water  
proof**

100% waterproof


**auto  
off**

Automatic switch off


**including  
protective  
cover**

With protective cover



Digital thermometer

Large Display

3-digit display

Audible signal when measuring stops

1 memory place

Measurement accuracy  $\pm 0.1^{\circ}\text{C}$  (between  $35.5$  and  $42^{\circ}\text{C}$ )Measurement accuracy  $\pm 0.2^{\circ}\text{C}$  (between  $32 - 35.4^{\circ}\text{C}$  and  $42.1 - 42.9^{\circ}\text{C}$ )

100% waterproof

Can be disinfected

Without glass, mercury free

Incl. 1 x 1,5V battery LR41

Battery changing indicator

Product measurements: 129 x 20 x 17 mm

Product weight: approx. 10 g ( including battery)

Medical device

5 year guarantee

Sales unit: 12 / Shipping carton: 12 x 12

EAN no.: 4211125 79116 2

Item no.: 791.16