



FT 85 Non-contact thermometer


 german
 engineering


Contactless infrared clinical thermometer
 Measures without skin contact - hygienic and safe



Laughing smiley - everything is OK /
 sad smiley - $\geq 38^{\circ}\text{C}$ temperature alarm



Measurement of forehead, object and room temperature in seconds

Measurement range:

Forehead temperature measurement: 34°C to $42,2^{\circ}\text{C}$ ($93,2^{\circ}\text{F}$ to 108°F)

Object temperature measurement: 0°C to 80°C (32°F to 176°F)

Measurement accuracy:

Forehead temperature measurement: $\pm 0,2^{\circ}\text{C}$ ($\pm 0,4^{\circ}\text{F}$)

Object temperature measurement: $\pm 4\%$ oder $\pm 2^{\circ}\text{C}$ ($\pm 4^{\circ}\text{F}$)

Measuring distance 2 - 3 cm

60 memory spaces (forehead measuring)

$^{\circ}\text{C}/^{\circ}\text{F}$ can be set

Display illumination

Sound control (On/Off)

Automatic switch off

Date and time

With battery changing indicator

Incl. 2 x 1.5V batteries AAA

Product measurements (LxWxH): ca. 161 x 42.8 x 36 mm

Product weight: approx. 67 g

Medical device

5 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 4 211125 79506 1

Item no.: 795.06