



measured without skin contact – hygiene and care



Laughing smiley – everything is OK /
 sad smiley – $\geq 38^{\circ}\text{C}$ temperature alarm



Measurement of forehead,
 object and room temperature in
 seconds



Particularly suitable for taking the temperature of babies and small children

60 memory spaces

$^{\circ}\text{C}/^{\circ}\text{F}$ can be set

Individually adjustable display illumination

Sound control (On/Off)

Automatic switch off

Includes storage box

With battery changing indicator

2 x 1,5V AAA batteries included

Product measurements: 188 x 47.6 x 29 mm

Product weight: approx. 82 g

Medical device

5 year guarantee

Sales unit: 4 / Shipping carton: 8x4

EAN no.: 4211125 79531 3

Item no.: 795.31