

living



germanengineering

HK 58 Cosy Rücken-Nacken Heizkissen

Flauschig weicher Heizkörper aus Micro-Flauschfaser

Atmungsaktiv, anschmiegsam und hautsympatisch

Körpergerechte Form für Rücken und Nacken

Innovativer LED Tastenschalter

Abnehmbare Schalter

Maschinenwaschbar bis 30°C

BSS (Beurer Sicherheits System)

Elektronische Temperaturregelung

6 beleuchtete Temperaturstufen

Schnellheizung

Abschaltautomatik nach ca. 90 Minuten

Abnehmbare Micro-Flauschfaser Bezug - waschbar

Öko-Tex 100

Die bei diesem Gerät eingesetzten Textilien erfüllen die hohen humanökologischen Anforderungen des Öko-Tex Standards 100, wie durch das Forschungsinstitut Hohenstein nachgewiesen.

Produktmaße: 62 x 42 cm

Produktgewicht: ca. 580 g

100 Watt

3 Jahre Garantie

VE: 5 / Transportkarton:

EAN-Nr.: 4211125 21505 7

Art.-Nr.: 215.05



HK 58 Cosy Back and Neck Heating pad

Cosy and soft heating pad made of cosy micro fibre

Breathable, cosy and skinfriendly

Formfitted form for back and neck

Innovative button switch with LED lights

Detachable control unit

Machine washable up to 30 °C

BSS (Beurer Safety System)

Electronic temperature control

6 illuminated temperature settings

Rapid heat up

Automatic switch off after approx. 90 minutes

Removable cosy micro fibre cover - washable

Oeko-Tex 100

The textiles used for this device meet the stringent human ecological requirements of Oeko-Tex Standard 100, as verified by Hohenstein Research Institute.

Product measurements: 62 x 42 cm

Product weight: approx. 580 g

100 watt

3 year guarantee

Sales unit: 5 / Shipping carton:

EAN no.: 4211125 21505 7

Item no.: 215.05



beurer