



BM 27+

beurer
medical

BM 27+ Blood pressure monitor



With USB-C connection for optional mains operation



Cuff position control*:
For ensuring the cuff is correctly
positioned on the upper arm



Scan & Save – the simple way
to save blood pressure
values in the app



Fully automatic blood pressure and
pulse measurement on the upper arm

Clear, easy-to-read display

Risk indicator: classification of measurements with coloured
scale

Arrhythmia detection: Warning in case of possible heart
rhythm disturbance

Four user memories for 30 measured values each

Average of all measurements

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. 4x1,5 V batteries AA and storage bag

Mains adapter available separately: Item no.: 072.78

Mains operation: Input: 100-240 V; 50-60 Hz; max. 0.5 A
Output: 5V DC, 1A

Medical device according to new Medical Device Regulation
(EU) 2017/745 (MDR)

Product measurements (LxWxH): 112 x 110 x 53 mm

Product weight: approx. 344.3 g (with cuff, without batteries)

5 year guarantee

Sales unit: 4 / Shipping carton: 6 x 4

EAN no.: 4211125 60418 9

Item no.: 604.18

*Cuff wrapping guide



604.18_BM27+_2024-06-21_03_DS_EN_BEU Subject to errors and changes